

Product Datasheet

Mouse IL13RA2 ELISA Kit (orb2991742)

Catalog Number	orb2991742
Category	Assays and Kits
Description	Mouse IL13RA2 ELISA Kit (96 Tests). Quantitate Mouse IL13RA2 in cell culture supernatants, serum, plasma(heparin, EDTA). Sensitivity: 10pg/ml. The brand indicates this is a premium quality ELISA kit. Each kit delivers precise quantification, high sensitivity, and excellent reproducibility. Only our most reliable and effective kits qualify as , guaranteeing top-tier results for your assays.
Reactivity	Mouse
UniProt ID	O88786
MW	44483 Da
Tested applications	ELISA
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	39.0-2,500 pg/mL
Sensitivity	10 pg/mL
Sample Types	cell culture supernatants, serum, plasma(heparin, EDTA)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)