

## Product Datasheet

### Mouse PRF1 (Perforin 1) ELISA Kit (orb2992016)

|                        |  |
|------------------------|--|
| <b>Catalog Number</b>  | orb2992016   |
| <b>Category</b>        | Assays and Kits  |
| <b>Description</b>     | Mouse PRF1 (Perforin 1) ELISA Kit  |
| <b>Target</b>          | PRF1   |
| <b>Reactivity</b>      | Mouse  |
| <b>UniProt ID</b>      | <b>P10820</b>  |
| <b>Assay Type</b>      | Sandwich   |
| <b>Assay Time</b>      | 120 minutes  |
| <b>Range</b>           | 0.156-10ng/ml  |
| <b>Sensitivity</b>     | 0.094ng/ml   |
| <b>Sample Types</b>    | Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples |
| <b>Storage</b>         | 2 - 8 °C for 12 months   |
| <b>Note</b>            | For research use only  |
| <b>Expiration Date</b> | 6 months from date of receipt.   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)