

Product Datasheet

Recombinant Human EDIL3 Protein (orb2992995)

Catalog Number	orb2992995
Category	Proteins
Description	Recombinant Human EGF-Like Repeat and Discoidin I-Like Domain-Containing Protein 3 is produced by our Mammalian expression system and the target gene encoding Asp24-Glu480 is expressed with a human IgG1 Fc tag at the C-terminus.
Conjugation	Unconjugated
Tag	C-Fc
Form/Appearance	Lyophilized
Buffer/Preservatives	0.2 µm filtered solution of 20mM Tris-HCl, 300mM NaCl, pH 8.0.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Purity	SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.
UniProt ID	O43854
MW	Predicted: 78.3 kDa. Observed: 79-97 kDa, reducing conditions
Expression System	Human cells
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Asp24-Glu480

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)