

Product Datasheet

Biotinylated Human ACE-2 (orb2993139)

Catalog Number	orb2993139
Category	Proteins
Description	Biotinylated Recombinant Human Angiotensin-Converting Enzyme 2 is produced by our Mammalian expression system and the target gene encoding Gln18-Ser740 is expressed with a Avi, 6His tag at the C-terminus.
Conjugation	Biotin
Tag	C-Avi-6His
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 300mM NaCl, 1mM ZnCl ₂ , 10% Glycerol, pH 7.4.
Purity	SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.
UniProt ID	Q9BYF1
MW	Predicted: 86.2 KDa. Observed: 110-140 KDa, reducing conditions
Expression System	Human Cells
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Gln18-Ser740
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)