

Product Datasheet

Human CASP10 (orb2993345)

Catalog Number	orb2993345
Category	Proteins
Description	Recombinant Human Caspase-10 is produced by our E.coli expression system and the target gene encoding Val220-Ile480 is expressed with a 6His tag at the C-terminus.
Conjugation	Unconjugated
Tag	C-6His
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 8% Sucrose, 1mM DTT, 0.05% Tween80, pH 8.5.
Purity	SDS-PAGE: Greater than 70% as determined by reducing SDS-PAGE. (QC verified)
UniProt ID	Q92851
MW	Predicted: 30.1 KDa. Observed: 33 KDa, reducing conditions
Expression System	E.coli
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Val220-Ile480
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)