

Product Datasheet

Human UBE2J2 (orb2993401)

Catalog Number	orb2993401
Category	Proteins
Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 J2 is produced by our E.coli expression system and the target gene encoding Glu124-Arg224 is expressed with a GST tag at the N-terminus.
Conjugation	Unconjugated
Tag	N-GST
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
Purity	SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.
UniProt ID	Q8N2K1
MW	Predicted: 37.96 KDa. Observed: 30&36 KDa, reducing conditions
Expression System	E.coli
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Glu124-Arg224
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558