

Product Datasheet

Human ATG3 (orb2993402)

Catalog Number	orb2993402
Category	Proteins
Description	Recombinant Human Ubiquitin-like-conjugating Enzyme ATG3/ATG3 is produced by our E.coli expression system and the target gene encoding Met1-Met314 is expressed.
Conjugation	Unconjugated
Tag	Tag free
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Purity	SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.
UniProt ID	Q9NT62
MW	Predicted: 36 KDa. Observed: 38 KDa, reducing conditions
Expression System	E.coli
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Met1-Met314
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)