

Product Datasheet

Recombinant Human EphA1 Protein (orb2994222)

Catalog Number	orb2994222
Category	Proteins
Description	Recombinant Human Ephrin A Receptor 1 is produced by our Mammalian expression system and the target gene encoding Lys26-Glu547 is expressed with a 6His tag at the C-terminus.
Conjugation	Unconjugated
Tag	C-6His
Form/Appearance	Lyophilized
Buffer/Preservatives	0.2 µm filtered solution of PBS, pH 7.4.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Purity	SDS-PAGE: Greater than 95% as determined by reducing SDS-PAGE.
UniProt ID	AAI30292.1
MW	Predicted: 57.4 KDa. Observed: 70-85 KDa, reducing conditions
Expression System	Human Cells
Endotoxins	1 EU/µg as determined by LAL test.
Biological Origin	Human
Expression Region	Lys26-Glu547

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.

Note

For research use only

Expiration Date

6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)