

Product Datasheet

Human IL-12 p40 ELISA Kit (orb2995300)

Catalog Number	orb2995300
Category	Assays and Kits
Description	Human IL-12 p40 ELISA Kit
Target	Human IL-12 p40
Reactivity	Human
UniProt ID	P29460
Tested applications	ELISA
Assay Type	Sandwich
Range	12.29 pg/ml - 3000 pg/ml
Sample Types	serum_plasma_cell_culture_supernatant
Storage	Shelf life: 6 months from date of receipt. Storage: 2 - 8 °C
Note	For research use only
NCBI	NP_002178.2
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)