

## Product Datasheet

### Human IL-5 ELISA Kit (orb2995301)

<b>Catalog Number</b>	orb2995301
<b>Category</b>	Assays and Kits
<b>Description</b>	Human IL-5 ELISA Kit
<b>Target</b>	Human IL-5
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>P05113</b>
<b>Tested applications</b>	ELISA
<b>Assay Type</b>	Sandwich
<b>Range</b>	2.74 pg/ml - 2000 pg/ml
<b>Sample Types</b>	serum_plasma_cell_culture_supernatant
<b>Storage</b>	Shelf life: 6 months from date of receipt. Storage: 2 - 8 °C
<b>Note</b>	For research use only
<b>NCBI</b>	<b>NP_000870.1</b>
<b>Expiration Date</b>	6 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)