

## **Product Datasheet**

## Human EDIL3(EGF Like Repeats And Discoidin I Like Domains Protein 3) Microsample ELISA Kit (orb2997687)

Catalog Number orb2997687

Category Assays and Kits

**Description** Human EDIL3(EGF Like Repeats And Discoidin I Like Domains Protein 3)

Microsample ELISA Kit

UniProt ID 043854

**Application notes** The test principle applied in this kit is Sandwich enzyme immunoassay. The

microtiter plate provided in this kit has been pre-coated with an antibody specific to Human EDIL3. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Human EDIL3. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Human EDIL3, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$  10nm. The

concentration of Human EDIL3 in the samples is then determined by comparing

the OD of the samples to the standard curve.

Assay Type Sandwich

Assay Time 3.5h

**Range** 0.32-20 ng/mL

**Sensitivity** 0.12 ng/mL

Sample Types Tissue homogenates, cell lysates and other biological fluids

For research use only

Note byt Ltd.

7 Signet Court, Swann Road Cambridge CB5 8LA

United Kingdom
Email: <a href="mailto:info@biorbyt.com">info@biorbyt.com</a>, <a href="mailto:support@biorbyt.com">support@biorbyt.com</a>

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC** 

68 TW Alexander Drive Research Triangle Park

Durham NC 27713 United States

Email: <u>info@biorbyt.com</u>, <u>support@biorbyt.com</u> Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558