

Product Datasheet

Goat Anti-Rat IgG H&L, BF405 conjugated (orb2999650)

Catalog Number	orb2999650
Category	Antibodies
Description	Goat Anti-Rat IgG H&L, BF405 conjugated
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	BF405
Reactivity	Rat
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein G
Immunogen	Native Rat IgG full length
Tested applications	FC, ICC, IF
Dilution range	ICC/IF=1:100-1000, IF=1:100-1000, Flow-Cyt=1:100-1000
Antibody Type	Secondary Antibody
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)