

Product Datasheet

Recombinant human IL5 protein, N-His (orb2999691)

Catalog Number	orb2999691
Category	Proteins
Description	Recombinant human IL5 protein, N-His
Target	Interleukin-5
Tag	N-His
Concentration	>0.5mg/ml
Purity	>90% as determined by SDS-PAGE
Purification	AC
Protein Sequence	20-134/134
UniProt ID	P05113
MW	15.5 kDa
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Expression System	E.coli
Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)