

Product Datasheet

IFIH1 Recombinant Rabbit mAb, FITC conjugated (orb3001189)

Catalog Number	orb3001189
Category	Antibodies
Description	IFIH1 Recombinant Rabbit mAb, FITC conjugated is an FITC conjugated antibody targeting IFIH1. This antibody is suitable for ICC. It exhibits reactivity with Human samples.
Target	IFIH1
Clonality	Recombinant
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Human
Buffer/Preservatives	10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% rAlbumin, 0.02% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein A
Immunogen	A synthesized peptide derived from human MDA5 (950-1025aa)
Tested applications	ICC
Dilution range	ICC/IF=1:50-200
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	5B1
Storage	Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)