

Product Datasheet

Human CCL17/TARC ELISA Kit (orb3001831)

Catalog Number	orb3001831
Category	Assays and Kits
Description	Human CCL17/TARC ELISA Kit
Reactivity	Human
Sensitivity	15.6-1000pg/ml
Sample Types	Cell culture supernates, serum, and plasma.
Storage	Store at 4°C for 12 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)