

Product Datasheet

Mouse Proopiomelanocortin (POMC) ELISA KIT (orb3012403)

Catalog Number	orb3012403
Category	Assays and Kits
Description	Mouse Proopiomelanocortin (POMC) ELISA KIT
Assay Type	Competitive
Range	0.15-10ng/mL
Sensitivity	0.07ng/mL
Sample Types	Serum;plasma and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)