

Product Datasheet

Rat Transient Receptor Potential Cation Channel Subfamily V, Member 1 (TRPV1) ELISA Kit (orb3012693)

Catalog Number	orb3012693
Category	Assays and Kits
Description	Rat Transient Receptor Potential Cation Channel Subfamily V, Member 1 (TRPV1) ELISA Kit
Reactivity	Rat
UniProt ID	O35433
Assay Type	Sandwich
Range	78.1-5000pg/mL
Sensitivity	29pg/mL
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)