

Product Datasheet

SAP97 Antibody: FITC (orb3013007)

Catalog Number	orb3013007
Category	Antibodies
Description	Rabbit Anti-Rat SAP97 Recombinant Monoclonal IgG Kappa
Clonality	Recombinant
Species/Host	Rabbit
Isotype	IgG Kappa
Conjugation	FITC
Reactivity	Human, Mouse, Rat
Concentration	1mg/mL
Buffer/Preservatives	640.91mM DMSO, 136.36 mM Ethanolamine, 126.89 mM chlorides, 9.09mM phosphates, 9.09mM NaHCO ₃
Purification	MabSelect Prisma
Immunogen	Fusion protein, AA #1-104 (N-terminus) of rat SAP97.
UniProt ID	Q62696
Tested applications	ICC, WB
Dilution range	WB (1:1000), ICC/IF (1:100); optimal dilutions for assays should be determined by the user.; optimal dilutions for assays should be determined by the user.
Application notes	1 µg/ml was sufficient for detection of SAP97 on Rat Brain Membrane lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Specificity	Detects N-terminal region of SAP97. Does not cross-react against PSD95, Chapsyn-110, or SAP102
Antibody Type	Recombinant Antibody
Clone Number	G28
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)