

Product Datasheet

Cyclooxygenase 1 Rabbit mAb Antibody (orb3015588)

Catalog Number	orb3015588
Category	Antibodies
Description	Cyclooxygenase 1 Rabbit mAb Antibody
Clonality	Monoclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
UniProt ID	P23219
Tested applications	ICC, IHC-Fr, IHC-P, IP, WB
Dilution range	WB: 1/500-1/1000, IHC: 1/50-1/100, IF: 1/50-1/200, IP: 1/20
Clone Number	R02-6G4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)