

Product Datasheet

Cyclosporin A-BSA (orb3021637)

Catalog Number	orb3021637
Category	Small Molecules
Description	Cyclosporin A-BSA is a hapten-carrier protein conjugate formed by coupling Cyclosporin A to bovine serum albumin (BSA). As a hapten, Cyclosporin A alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. Cyclosporin A-BSA may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01 M PBS (pH 7.4) buffer at a typical concentration of 1 mg/mL.
Target	Others
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)