

Product Datasheet

DL-Serine (Standard) (orb3025614)

Catalog Number	orb3025614
Category	Small Molecules
Description	DL-Serine (Standard) is the standard substance of DL-Serine, and it is applicable for quantitative analysis, quality control, and related research in biochemical experiments. DL-Serine is a non-essential amino acid and a natural ligand and allosteric activator of pyruvate kinase M2.
MW	105.09
CAS Number	302-84-1
Formula	C ₃ H ₇ NO ₃
SMILES	NC(CO)C(O)=O
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)