

## Product Datasheet

# Polyclonal Rabbit Desmin K5 Antibody (orb302935)

<b>Catalog Number</b>	orb302935
<b>Category</b>	Antibodies
<b>Description</b>	Rabbit polyclonal antibody to Desmin K5
<b>Clonality</b>	Polyclonal
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Reactivity</b>	Bovine, Gallus, Goat, Hamster, Human, Mouse, Rat, Zebrafish
<b>Buffer/Preservatives</b>	PBS with 0.09% sodium azide
<b>Immunogen</b>	Please contact us for more information.
<b>Tested applications</b>	IHC-Fr, IHC-P, WB
<b>Dilution range</b>	IHC-fr: 1:25-1:100, IHC-P: 1:25-1:100, WB: 1:100-1:500
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)