

Product Datasheet

Hedgehog/Smoothened,Apoptosis Antibody (orb3031945)

Catalog Number	orb3031945
Category	Antibodies
Description	ERY-974 is a T cell-redirecting bispecific antibody (TCB) targeting GPC3. It binds to both GPC3 and the CD3 receptor on the surface of T cells, activating them to release perforin and granzyme, which induces apoptosis in tumor cells. ERY-974 shows potential for research in solid tumors, such as liver cancer.
Target	Apoptosis,Hedgehog/Smoothened
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)