

Product Datasheet

PD-1/PD-L1, VEGFR Antibody (orb3031976)

Catalog Number	orb3031976
Category	Antibodies
Description	Pumitamig is a humanized bispecific immunoglobulin (H- γ 1-VH _L - κ) dimer antibody that targets human vascular endothelial growth factor A (VEGFA) and human programmed death-ligand 1 (PD-L1). It holds promise for research in various solid tumors.
Target	PD-1/PD-L1, VEGFR
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2999726-97-3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)