

Product Datasheet

Murine RNase Inhibitor, GMP-Grade (40 U/ μ L) (orb3032370)

Catalog Number	orb3032370
Category	Tools
Description	Murine RNase Inhibitor, GMP-Grade is a high-purity enzyme protecting RNA from degradation. Ideal for diagnostics, therapeutics, and research. GMP-compliant.
Storage	Store at -25°C ~ -15°C for up to 1 year.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)