

Product Datasheet

DNase I, GMP-Grade (2 U/ μ L) (orb3032374)

| | |
|------------------------|---|
| Catalog Number | orb3032374 |
| Category | Tools |
| Description | DNase I, GMP-Grade (2 U/ μ L) is a high-purity enzyme for DNA removal in RNA and protein samples, ensuring purity for RT-PCR, qPCR, and other applications. |
| Storage | Store at -25°C ~ -15°C for up to 1 year. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)