

## Product Datasheet

### DNase I, GMP-Grade (2 U/ $\mu$ L) (orb3032374)

<b>Catalog Number</b>	orb3032374
<b>Category</b>	Tools
<b>Description</b>	DNase I, GMP-Grade (2 U/ $\mu$ L) is a high-purity enzyme for DNA removal in RNA and protein samples, ensuring purity for RT-PCR, qPCR, and other applications.
<b>Storage</b>	Store at -25°C ~ -15°C for up to 1 year.
<b>Note</b>	For research use only

---

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)