

Product Datasheet

ssDNA Assay Kit (orb3032376)

Catalog Number	orb3032376
Category	Assays and Kits
Description	ssDNA Assay Kit provides fast, sensitive ssDNA quantification with 1-200 ng range, high contaminant tolerance, and compatibility with Qubit & microplate readers.
Storage	Ship with gel packs. Store at 2-8°C, protected from light, for up to 6 months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)