

Product Datasheet

RNA Cleaner (orb3032379)

Catalog Number	orb3032379
Category	Tools
Description	RNA Cleaner ensures efficient RNA purification with magnetic beads, removing proteins and impurities, which can be used for NGS, RT-PCR, qRT-PCR, and more.
Storage	Store at 2-8°C for up to 1 year. Avoid freezing!
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)