

Product Datasheet

Universal TaqMan qPCR Mix (UDG Plus) (orb3032389)

Catalog Number	orb3032389
Category	Tools
Description	Universal TaqMan qPCR Mix (UDG Plus) is a multiplex fluorescence-based qPCR premix with a hot start Taq enzyme and dUTP/UDG anti-contamination system.
Storage	Store at -25°C ~ -15°C for up to 1 year. Ship on dry ice.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)