

Product Datasheet

Universal Blue qPCR Master Mix (Dye-Based) (orb3032391)

Catalog Number	orb3032391
Category	Tools
Description	Universal Blue qPCR Master Mix (Dye-Based) is a premix for real-time quantitative PCR (qPCR) with hot start Taq DNA polymerase for highly sensitive reaction.
Storage	Store at -15°C to -25°C for up to 18 months. Ship with ice packs. This product contains fluorescent dyes, so it should be protected from strong light during storage and while preparing the reaction system.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)