

Product Datasheet

dNTP Mix (25 mM each) (orb3032410)

Catalog Number	orb3032410
Category	Tools
Description	dNTP Mix (dATP, dCTP, dGTP, dTTP) with 99% purity (HPLC) and concentration 25 mM for each dNTP, perfect for PCR, cDNA synthesis, DNA sequencing, and more.
Storage	Store at -25°C ~ -15°C for up to 24 months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)