

Product Datasheet

1st Strand cDNA Synthesis SuperMix for qPCR (gDNA Digester Plus) (orb3032414)

Catalog Number	orb3032414
Category	Tools
Description	1st Strand cDNA Synthesis SuperMix for qPCR (gDNA Digester Plus) with Reverse Transcriptase, High-Yield, Fast, Thermostable is designed for two-step RT-qPCR.
Storage	Store at -25°C ~ -15°C for up to 18 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)