

Product Datasheet

Human Complement Component 5a Receptor 1 (C5aR1) ELISA Kit (Ready to Use) (orb3032633)

| | |
|------------------------|--|
| Catalog Number | orb3032633 |
| Category | Assays and Kits |
| Description | This Human Complement Component 5a Receptor 1 (C5aR1) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in tissue homogenates, cell lysates and other biological fluids . This product is for Research Use Only. |
| Reactivity | Human |
| UniProt ID | P21730 |
| Assay Type | Sandwich |
| Assay Time | 3 hours |
| Range | 78.125-5000 pg/mL |
| Sensitivity | 36.5 pg/mL |
| Sample Types | tissue homogenates, cell lysates and other biological fluids |
| Storage | Store at -20°C |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Standard curve for orb3032633 ELISA kit.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)