

Product Datasheet

FS-Human PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit (orb3033230)

Catalog Number	orb3033230
Category	Assays and Kits
Description	FS-Human PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
Target	PECAM1/CD31
Reactivity	Human
UniProt ID	P16284
Assay Type	Sandwich
Assay Time	1.5H
Standard	10ng/mL
Range	0.16-10ng/mL
Sensitivity	0.03 ng/mL
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)