

Product Datasheet

ACVRL1 (NM_000020) Human Recombinant Protein (orb3033934)

Catalog Number	orb3033934
Category	Proteins
Description	Purified recombinant protein of Human activin A receptor type II-like 1 (ACVRL1), transcript variant 1, Asp22-Leu119, with N-terminal His-ABP tag, expressed in E. coli, 50ug
Tag	N-His-ABP
Buffer/Preservatives	25mM Tris, pH 8.0, 150 mM NaCl, 10% glycerol
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
UniProt ID	P37023, A0A0S2Z310
MW	25.9 kDa
Expression System	E. coli
Biological Origin	Human
Storage	Store at -80°C.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)