

## Product Datasheet

### PLXDC1 peptide (orb304750)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block Anti-PLXDC1 Antibody (orb304509) reactivity.                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                 |
| <b>Tested Applications</b> | Blocking                                                                                                                                                                          |
| <b>Target</b>              | PLXDC1                                                                                                                                                                            |
| <b>Form/Appearance</b>     | Lyophilized powder                                                                                                                                                                |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Application notes</b>   | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| <b>Purity</b>              | > 85%                                                                                                                                                                             |
| <b>Source</b>              | Synthetic                                                                                                                                                                         |
| <b>Uniprot ID</b>          | <b>Q8IUK5, Q91ZV7</b>                                                                                                                                                             |
| <b>Entrez</b>              | <b>57125, 72324</b>                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)