

Product Datasheet

Hyaluronidase PH20 (SPAM1) (NM_003117) Human Recombinant Protein (orb3049771)

Catalog Number	orb3049771
Category	Proteins
Description	Recombinant protein of human sperm adhesion molecule 1 (PH-20 hyaluronidase, zona pellucida binding) (SPAM1), transcript variant 1
Tag	C-Myc/DDK
Buffer/Preservatives	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
UniProt ID	P38567, Q5D1J4
MW	58.2 kDa
Expression System	HEK293T
Biological Origin	Human
Storage	Store at -80°C.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)