

Product Datasheet

Human Hyaluronidase PH-20(SPAM1) ELISA kit (orb3053039)

Catalog Number	orb3053039
Category	Assays and Kits
Description	Human Hyaluronidase PH-20(SPAM1) ELISA kit
Target	sperm adhesion molecule 1 (PH-20 hyaluronidase, zona pellucida binding)
Reactivity	Human
UniProt ID	P38567
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	78 pg/mL-5000 pg/mL
Sensitivity	39 pg/mL
Sample Types	cell culture supernates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)