

Product Datasheet

Monkeypox Virus (MPXV) A35R Humanized Monoclonal Antibody (Detection) (orb3054735)

Catalog Number	orb3054735
Category	Antibodies
Description	Monkeypox Virus (MPXV) A35R Humanized Monoclonal Antibody (Detection)
Clonality	Monoclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human Humanized
Form/Appearance	Purified,Liquid
Buffer/Preservatives	PBS, 5% mannitol, 0.2% Tween20, pH7.4
Purity	>95% pure (SDS-PAGE)
MW	150 kDa
Tested applications	CLIA
Specificity	Specific for MPXV A35R
Matched Pairs	Control:orb3054736; Capture: orb3054734; Detection: orb3054735;
Source	Mammalian cell
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)