

Product Datasheet

Recombinant Helicobacter Pylori (HP) KatA Protein (Control) (orb3054743)

| | |
|-----------------------------|--|
| Catalog Number | orb3054743 |
| Category | Proteins |
| Description | Recombinant Helicobacter Pylori (HP) KatA Protein (Control) |
| Form/Appearance | Purified, Liquid |
| Buffer/Preservatives | 50mM Tris-HCl, 500mM N/ACl, 5% Glycerol, pH 8.0 |
| Purity | >95% pure (SDS-PAGE) |
| MW | 58.63 kDa |
| Tested applications | CLIA |
| Expression System | E. coli |
| Matched Pairs | Control: orb3054743; Capture: orb2280039; Detection: orb2280038; |
| Source | Helicobacter Pyloritryptococcus |
| Storage | For long term storage, the product should be stored at -20°C or lower. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)