

Product Datasheet

SARS COV2 Nucleoprotein (NP) (His Tag) Protein (Control) (orb3054751)

Catalog Number	orb3054751
Category	Proteins
Description	SARS COV2 Nucleoprotein (NP) (His Tag) Protein (Control)
Form/Appearance	Purified,Liquid
Buffer/Preservatives	PBS, 5% mannital, 0.01% Tween 80, pH7.4
Purity	>95% pure (SDS-PAGE)
MW	47 kDa
Tested applications	CLIA
Expression System	Mammalian cell
Matched Pairs	Control: orb3054751; Capture: orb3054750; Detection: orb3054749;
Source	Severe acute respiratory syndrome coroN/Avirus 2 (2019-nCoV) (SARS-CoV-2)
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)