

Product Datasheet

Adenovirus (ADV) Mouse Monoclonal Antibody (Capture) (orb3054755)

Catalog Number	orb3054755
Category	Antibodies
Description	Adenovirus (ADV) Mouse Monoclonal Antibody (Capture)
Clonality	Monoclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Purified, Liquid
Buffer/Preservatives	PBS, pH 7.4
Purity	>90% pure (SDS-PAGE).
MW	150 kDa
Tested applications	CLIA
Specificity	Specific for Adenovirus
Matched Pairs	Control: orb1921798; Capture: orb3054755; Detection: orb1921796;
Source	Mammalian cell
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)