

## Product Datasheet

# Recombinant Human CD114/CSF3R Protein, C-His (orb3059442)

|                          |                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb3059442                                                                                                                                                                |
| <b>Category</b>          | Proteins                                                                                                                                                                  |
| <b>Description</b>       | Recombinant Human CD114/CSF3R Protein, C-His                                                                                                                              |
| <b>Tag</b>               | His-tag                                                                                                                                                                   |
| <b>Expression System</b> | Mammalian Cells                                                                                                                                                           |
| <b>Storage</b>           | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| <b>Note</b>              | For research use only                                                                                                                                                     |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                            |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)