

Product Datasheet

FGF2 (Detection) Mouse Monoclonal Antibody (orb3073320)

| | |
|-----------------------------|---|
| Catalog Number | orb3073320 |
| Category | Antibodies |
| Description | Mouse monoclonal detection antibody to FGF2 |
| Target | FGF2 |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Isotype | Mouse IgG1 |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Buffer/Preservatives | 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide. |
| Purification | The antibody was purified by Protein A/G affinity chromatography. |
| Immunogen | Recombinant protein of human FGF2 |
| UniProt ID | P09038 |
| Tested applications | ELISA |
| Specificity | Recognizes native and recombinant human FGF2 protein |
| Matched Pairs | Capture: orb3073321; Detection: orb3073320 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Entrez | 2247 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)