

Product Datasheet

Human Cys-C(Cystatin C) ELISA Kit (orb3074068)

Catalog Number	orb3074068
Category	Assays and Kits
Description	Human Cys-C(Cystatin C) ELISA Kit
Reactivity	Human
Assay Type	Sandwich
Assay Time	3.5h
Range	3.13-200 ng/mL
Sensitivity	2.47 ng/mL
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)