

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 1B protein (C-His) (orb312061)

Catalog Number	orb312061
Category	Proteins
Description	Recombinant Human Tumor Necrosis Factor Receptor II is produced by our Mammalian expression system and the target gene encoding Leu23-Asp257 is expressed with a 6His tag at the C-terminus.
Tag	C-His
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
UniProt ID	P20333
Application notes	6His tag at the C-terminus.
Expression System	Human Cells
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Biological Origin	Human
Expression Region	Leu23-Asp257
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558