

Product Datasheet

Mouse Chemokine C-X3-C-Motif Receptor 1 (CX3CR1) ELISA Kit (orb3124139)

| | |
|------------------------|--|
| Catalog Number | orb3124139 |
| Category | Assays and Kits |
| Description | Mouse Chemokine C-X3-C-Motif Receptor 1 (CX3CR1) ELISA Kit |
| Reactivity | Mouse |
| Range | 0.156-10ng/mL |
| Sensitivity | 0.047ng/mL |
| Sample Types | Tissue homogenates, cell lysates or other biological fluids. |
| Storage | Store at 4°C & -20°C for 12 months. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)