

Product Datasheet

Recombinant human IL-6 protein (Active, CHO) (orb3124196)

Catalog Number	orb3124196
Category	Proteins
Description	Recombinant human IL-6 protein (Active, CHO)
Target	Interleukin-6
Tag	Tag free
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized from a 0.22µm-filtered solution containing PBS, 5% Mannitol and 0.01% Tween 80, pH 7.4
Reconstitution	It is recommended to redissolve in sterile deionized water
Purity	≥ 95% as determined by SDS-PAGE.
Purification	AC
Protein Sequence	30-212/212
UniProt ID	P05231
MW	20.8 kDa
Expression System	CHO cell
Endotoxins	≤10 EU/mg
Biological Origin	Human

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 34°C under sterile conditions after reconstitution.

Note

For research use only

Expiration Date

6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)