

Product Datasheet

tdTomato Antibody (orb3127246)

Catalog Number	orb3127246
Category	Antibodies
Description	tdTomato goat polyclonal antibody, which has been validated for WB and IF using 293T cells overexpressing tdTomato.
Species/Host	Goat
Conjugation	Unconjugated
Reactivity	All
Form/Appearance	Liquid
Concentration	~1mg/ml
Buffer/Preservatives	PBS, 50% glycerol, 0.05% Proclin 300, 0.08%rAlbumin
Purification	Ammonium sulfate precipitation&Antigen affinity chromatography
Immunogen	Recombinant tdTomato protein
Tested applications	IF, IHC-Fr, WB
Dilution range	WB: 1:500-1000, IHC-F: 1:100-500, IF: 1:100-300
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot analysis of extracts from 293T cells



Immunofluorescence analysis of HEK293 cells using tdTomato antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)